



LIFTONIC PLUS ONDA DINAMICA®

SELENIA ITALIA

OUR HISTORY

SELENIA ITALIA SRL is a Beauty Company founded in 1991 by the Posarelli family, a leader in the field of professional aesthetics and aesthetic medicine. The company is specialized in the production and commercialization of innovative technologies and hi-tech cosmetics and is well known for its reliability, seriousness and maximum support to its customers. Selenia Italia is a certified company that distributes its products in Italy and in many European, African and Middle Eastern countries.



MADE IN ITALY

Innovation, effectiveness, reliability and Made in Italy design distinguish all the production of Selenia Italia, which has always been committed to the continuous pursuit of excellence. At the headquarters in Pisa, a specialized staff, in collaboration with a scientific team made up of doctors, engineers and technicians, devote themselves every day to the study and creation of new technologies and products to offer customers real, tangible and persistent results.

LIFTONIC PLUS ONDA DINAMICA®

PURE REGENERATING ENERGY

Liftonic Plus Onda Dinamica® is the new technology designed to reactivate the vitality of the tissues and obtain immediate and persistent effects of regeneration and toning. Specific for firming, skin remodeling, evident reduction of wrinkles and facial imperfections.

ONDA DINAMICA® THE NEW PATENTED AND MULTI-ACTION FORM OF ENERGY

An exclusive and patented technology born from the experience and technical scientific research of Selenia Italia. This form of energy stimulates the tissue causing a cellular resonance that acts on the dystrophy and on the causes of imperfections thanks to the stimulation of slowed biochemical reactions and a prompt metabolic reactivation.





This innovative combination of energies represents the maximum expression of advanced technologies for the reduction of facial imperfections such as widespread relaxation, wrinkles, dull skin and deep atony.

SUPERLATIVE RESULTS

The "dynamic" delivery of the energies emitted by each handpiece allows for a powerful multiple stimulation that achieves, in a short time, effects of reduction, remodeling and firming of the tissues that are evident and persistent over time.

AN EFFECTIVE AND PLEASANT TREATMENT

With each application, the tissue is heated, vascularized, oxygenated, stimulated, defibrotized and sculpted simultaneously, thanks to the action of 4 combined energies, dosed to obtain a complete dermal-hypodermic restructuring, visible from the first session.

AN EXCLUSIVE, PATENTED AND CLINICALLY TESTED TECHNIQUE

Thanks to the Onda Dinamica® treatment it is possible to obtain, in each session, superior results compared to those obtained with other types of equipment, halving the metabolic activation time.

LIFTONIC PLUS ONDA DINAMICA®

ONE HANDPIECE FOR FACIAL ANALYSIS AND TWO SPECIFIC HANDPIECES TO AWAKEN THE SKIN'S VITALITY AND FIGHT AGAINST IT. CAUSES OF SKIN AGING.



SKIN ANALYSIS

It identifies skin alterations and determines the type of treatment to be carried out.



ONDA DINAMICA® FACE LIFT HANDPIECE

It has a marked anti-aging, vascularizing, decontracting, regenerating, defibrotizing and tonic / firming action.



DERMOTONIC FACE HANDPIECE

It allows the penetration and diffusion of cosmetics and active ingredients and stimulates osmotic exchanges and oxygenation, restoring elasticity and tone to the tissues.



A UNIQUE EQUIPMENT FOR TREATMENTS ON ANY SKIN TYPE:

- DEHYDRATED
- DISCHROMIC
- MATURE
- OILY / IMPURE
- SENSITIVE
- AGED
- RELAXED / DULL
- ACNEIC

FOR EVERY KIND OF FACE IMPERFECTION:

- SKIN AGING
- ATONY
- WRINKLES
- WATER RETENTION
- BAGS
- DARK CIRCLES



OndaDINAMICA®
SELÈNIA ITALIA


ONDA DINAMICA® FACE HANDPIECE LIFT

ANTI-AGING ACTION
DERMAL - HYPODERMIC

The Onda Dinamica® delivered by the regenerating handpiece allows the application of 4 energies at the same time and therefore to obtain a trophic and vascularizing effect: it activates the skin metabolisms and stimulates the endogenous production of collagen, elastin and hyaluronic acid, with a lifting, plumping and toning / firming action on the skin.

A powerful regenerative multiple stimulation that combines 4 energies: low frequency radiofrequency, high frequency ultrasound, infrared and red photostimulation delivered simultaneously. It defibrotizes the tissues, decreases inflammation, stimulates cell metabolisms and activates a quick recovery of slowed vitality, with evident results from the first session.

The skin is immediately oxygenated, elasticized, plumped, revitalized and toned. This treatment is suitable for all skin types and facial imperfections, even in combination with specific cosmetics.



IT ACTS ON:

1 EPIDERMIS

It stimulates the turnover of the surface layers with a lifting and anti-aging effect.

2 DERMIS

It has a regenerative action, stimulating the production of collagen. It defibrotizes, elasticizes, tones and remodels simultaneously.

3 SUBCUTANEOUS TISSUE

It increases the veno-lymphatic function, oxygenation and subcutaneous metabolisms.

4 MUSCULATURE

It increases the arteriovenous circulation, the relaxation of muscle fibers and the cellular metabolisms.




DERMOTONIC FACE HANDPIECE

DERMAL - EPIDERMAL ACTION FOR
THE DEFIBROTIZATION OF THE SKIN
AND THE VEHICULATION OF ACTIVE
INGREDIENTS

Dermoporation is an exclusive skin stimulation technique that delivers 4 energies simultaneously: it allows the penetration and diffusion of cosmetics and active ingredients into the skin, it stimulates osmotic exchanges and oxygenation, and restores elasticity and tone to the tissues.

It uses an innovative combination of micro currents with a carrier and tonic action in combination with ultrasonic sonophoresis and infrared and green photostimulation. An innovative treatment that provides the skin with a quantity of stimulating energies that allow it an immediate functional response and make it visibly toned, firm, luminous and hydrated.

This treatment is suitable for all skin types and facial imperfections, in combination with specific cosmetics.



1 IT VEHICULATES **WATERSOLUBLE ACTIVE INGREDIENTS** and increases their effectiveness.

2 IT VEHICULATES **LIPOSOLUBLE ACTIVE INGREDIENTS** and makes the skin more toned and plumped up.

3 IT **DEFIBROTIZES AND ELASTICIZES** the skin, favoring greater oxygenation.

4 IT **STIMULATES THE DRAINAGE AND ACTS ON SKIN INFLAMMATION**, reducing toxins and aging factors.



SKIN ANALYSIS

ANALYZE SKIN AND CALCULATE
BASED ON SPECIFIC ALGORITHMS
THE TYPE OF TREATMENT NECESSARY

The Analyzer has been designed to highlight the aesthetic needs of the customer and to give value to the operator: it simplifies the work and allows to be more effective in proposing treatments and self-care products. Also useful for monitoring the obtained results.

The result of the analysis is displayed on screen, printed on paper, or can be sent to the customer via e-mail.



How does it work?

The analysis of the skin takes place through a particular multifrequency bio-impedancemetry system that allows the determination of the amount of water and its distribution in the various intra and extra cellular districts.

Through this analysis, a visual examination of the surface imperfections and a brief investigation of the customer's habits, the device allows to trace a precise profile of the skin based on 15 parameters, such as: Total, intracellular and extracellular water, water retention, superficial dehydration, stratum corneum, sebum, tone, elasticity, sensitivity, spots, wrinkles and Ph, allowing the determination of the type and number of treatments to be performed, the beauty routine products and practical tips for maintaining healthy and beautiful skin.





THE ADVANTAGES OF THE DEVICE

BIO-FEEDBACK SAFETY SYSTEM

Bio-Feedback is the safety system that automatically monitors and adapts energy emissions basing on the biological response and the absorption capacity of the tissues. It allows to carry out effective treatments in absolute relaxation, without any harmful effects.

1.

TRAINING FOR USE

Selenia Italia provides training with experts and professionals of the sector as well as constantly updated training material and multimedia supports.

2.

ITALIAN DESIGN

The device was designed with the ingenuity and creativity that characterize Italian products. Its modern, practical, functional, and compact design fits elegantly into any beauty salon.

3.

10.1 " DISPLAY WITH HIGH RESOLUTION TOUCH SCREEN

4.

A practical and intuitive interface provides all the necessary information for personalized programming.

5.

MEN - WOMEN PROGRAMS AND LANGUAGE CHOICE

The protocols are designed for the use of the device on men and women, also in combination with other treatments. The software is available in: Italian, French, English, Spanish and German.

INTEGRATED SELF-DIAGNOSTIC

A specific control system verifies the efficiency of the device and the handpieces in real time, allowing to work in absolute safety.

6.

REAL-TIME ASSISTANCE

Thanks to the "wired Internet" technology, each device is constantly interconnected to the company, allowing remote assistance to be guaranteed in real time.

SELENIA NETWORK

Liftonic Plus Onda Dinamica® is a device with advanced technology, designed to be connected to the internet and external PC via WI-FI or LAN and can perform periodic upgrades, offering an always up-to-date device.

This application makes the company always connectable and available to its customers anywhere in the world:

- Software updates
- Equipment monitoring
- Remote technical assistance
- Archiving of data
- Download of new programs

7.

9.

SELENIA COMMUNITY

Becoming members of the Selenia Italia community, the users have the possibility to access the network and view directly on the monitor:

- Training video tutorials
- Online treatment protocols
- Updated procedures
- Products to use
- Advice and tips on treatments
- Personal reminders about treatments



OUR RESULTS





SAFETY AND QUALITY



EFFECTIVENESS GUARANTEE

Certified device for aesthetic use
In compliance with interministerial
decree n°205



QUALITY GUARANTEE

Selenia Italia is a company certified according
to the EN ISO 13485-2016 Quality System.



SOFTWARE

Display with high resolution graphic interface.
Proprietary software remotely upgradeable via Wi-Fi or via USB or RS-232 port.
Remote control and remote diagnostics.



DATA SHEET

Commercial name	LIFTONIC PLUS ONDA DINAMICA®
General power supply	110/230 Vac - 50/60 Hz
MAX absorbed power	200 VA
Radiofrequency Base frequency	500 KHz
Ultrasound 1 Base frequency	3400 KHz ± 10%
Electroporation Pulse frequency	200 Hz
IR LEDS	880 nm
GREEN LED	515 nm
RED LED	635 nm
Weight	15Kg
Dimensions (with trolley)	134x45x32 cm

Device Classification: type BF class I

Selenia Italia reserves the right to make changes to the equipment at any time and without notice. Selenia Italia equipment is designed and programmed to have results on imperfections without exceeding the permitted aesthetic limit.

Scientific studies

1. Manuel Albornoz-Cabello, Alfonso Javier Ibáñez-Vera & Blanca De la Cruz-Torres (2017) Efficacy of monopolar dielectric transmission radio frequency in panniculus adiposus and cellulite reduction, Journal of Cosmetic and Laser Therapy, 19:7, 422-426, DOI: 10.1080/14764172.2017.1342041
2. Ibanez-Vera AJ, Radiofrequency in aesthetics skin treatment: Classification and modalities, Journal dermatology research skin care, 2017, 1(1):1-3
3. R. Stephen Mulholland, Malcolm D. Paul, Charbel Chalfoun, Noninvasive Body Contouring with Radiofrequency, Ultrasound, Cryolipolysis, and Low-Level Laser Therapy, Clin Plastic Surg 38 (2011) 503–520, doi:10.1016 / j.cps.2011.05.002
4. Deepal V Wakade, Chitra S Nayak, Kalpana D Bhatt, A study comparing the efficacy of monopolar radiofrequency and glycolic acid peels in facial rejuvenation of aging skin using histopathology and ultrabiomicroscopic sonography (UBM) – An evidence based study, ACTA MEDICA (Hradec Králové) 2016; 59(1):14–17
5. El-Domyati M, el-Ammawi TS, Medhat W, et al. Radiofrequency facial rejuvenation: Evidence-based effect. J Am Acad Dermatol 2011 Mar; 64(3): 524–35
6. Sadick NS, Mulholland RS. A prospective clinical study to evaluate the efficacy and safety of cellulite treatment using the combination of optical and RF energies for subcutaneous tissue heating. J Cosmet Laser Ther 2004;6:187–90
7. Brightman L, Weiss E, Chapas AM, et al. Improvement in arm and post partum abdominal and flank subcutaneous fat deposits and skin laxity using a bipolar radiofrequency, infrared, vacuum and mechanical massage device. Lasers Surg Med 2009;41:791–8.
8. Anolik R, Chapas AM, Brightman LA, et al. Radiofrequency devices for body shaping: a review and study of 12 patients. Semin Cutan Med Surg 2009; 28:236–43.
9. Rubbani S. Advances in monopolar radiofrequency for the treatment of stretch marks in the arms, thighs and abdomen [abstracts]. American Society for Laser Medicine and Surgery 2008;370:111.
10. Rohrich RJ, Beran SJ, Kenkel JM, et al. Extending the role of liposuction in body contouring with ultrasound-assisted liposuction. Plast Reconstr Surg 1998; 101(4):1090–102 [discussion: 1117–9].
11. Teitelbaum SA, Burns JL, Kubota J, et al. Noninvasive body contouring by focused ultrasound: safety and efficacy of the Contour I device in a multicenter, controlled clinical study. Plast Reconstr Surg 2007;120(3):779–89.
12. Brown SA, Greenbaum L, Shtukmaster S, et al. Characterization of nonthermal focused ultrasound for non-invasive selective fat cell disruption (lysis): technical and preclinical assessment. Plast Reconstr Surg 2009;24(1):92–101.
13. Moreno-Moraga J, Valero-Altes T, Riquelme AM, et al. Body contouring by non-invasive transdermal focused ultrasound. Lasers Surg Med 2007;39:315–23
14. Ascher B. Safety and efficacy of UltraShape contour 1 treatments to improve the appearance of body contours: multiples treatments in shorter intervals. Aesthet Surg J 2010;30(2):217–24.
15. Leal H. Combined modality of focused ultrasound and radiofrequency for non-invasive fat disruption and body contouring—results of a single treatment session. Presented at IMCAS Paris, January 9, 2010.
16. Fatemi A, Kane MAC. High-intensity focused ultrasound effectively reduces waist circumference by ablating adipose tissue from the abdomen and flanks: a retrospective case series. Aesthetic Plast Surg 2010;34(5):577–82.
17. Kulick M. Evaluation of the combination of radio frequency, infrared energy and mechanical rollers with suction to improve skin surface irregularities (cellulite) in a limited treated area. J Cosmet Laser Ther 2006;8:185–90.
18. Neira R, Arroyave J, Ramirez H, et al. Fat liquefaction: effect of low-level laser energy on adipose tissue. Plast Reconstr Surg 2002;110(3):912–22.
19. Shinichiro Takezaki, Tokuya Omi, Shigeru Sato, and Seiji Kawana, Ultrastructural observation of human skin following irradiation with visible red light-emitting diodes (leds): a preliminary in vivo report, Laser Therapy 14.4: 153 – 160, Accepted for publication: December 2005
20. Karu TI. Low power laser therapy. Biomedical Photonics Handbooks 2003; CRC Press. C 48
21. Weiss RA et al. Clinical trial of a novel non thermal Led array for reversal of photoaging: clinical, histologic and surface profilometric results. Laser Surg Med 2005;36:85-91.
22. Karu T et al. A novel mitochondrial signalling pathway activated by visible – to near infrared radiation. Photochem Photobiol 2004;80(2): 366-372.

LIFTONIC PLUS
ONDA DINAMICA®



